

BLICKMAN



Surgical Utensils

Safe, sterile, stainless steel

Includes:

- Dressing Jars
- Emesis Basins
- Ointment Jars
- Bucket Pail
- Solution Basins
- Instrument Trays
- Catheter Tray
- Iodine Cups
- Medicine Cup
- Graduated Measuring Cups



THE BLICKMAN DIFFERENCE

- We offer a comprehensive range of surgical utensils to complement the high-end quality of our hospital equipment.
- We offer utensils that function to hold other surgical materials and liquids of various types.
- All utensils are stainless steel and welded smooth to ensure sterility for surgery.
- Our measuring utensils have units of measurement commonly used in the medical industry etched inside.

Key Product Features



Surgical Utensils

Description: Our surgical utensils are solid stainless steel and are great for sterile use in the OR.



Dressing Jars



Emesis Basins



Ointment Jars

Model #	Item #	Capacity/Dimensions	Shipping Weight
Dressing Jars			
DJ-1	9788000001	0.5 liters	1 lbs.
DJ-2	9788010002	1.0 liters	1 lbs.
DJ-3	9788020003	1.9 liters	2 lbs.
DJ-4	9788030004	2.8 liters	2 lbs.
DJ-5	9788040005	4.0 liters	2 lbs.
DJ-7	9788080007	7.6 liters	3 lbs.
Emesis Basins			
EB-1	9717090400	12 oz.	1 lbs.
EB-2	9720595450	16 oz.	1 lbs.
EB-3	9725511555	28 oz.	1 lbs.
EB-4	9730514065	56 oz.	1 lbs.
Ointment Jar			
OJ-1	9787720001	8 oz.	1 lbs.



Bucket Pail
KB-1



Solution Basins
SB-1

Model #	Item #	Capacity/Dimensions	Shipping Weight
Bucket Pail			
KB-1	9701920000	12.3 liters	4 lbs.
Solution Basin			
SB-1	9737370120	8 liters	3 lbs.
1830	9701830000	6.5 liters	2 lbs.

Surgical Utensils

Description: Our surgical utensils are solid stainless steel and are great for sterile use in the OR.



Instrument Trays
IT1015-2520



Catheter Tray
CT1, IT, IC

Model #	Item #	Capacity/Dimensions	Shipping Weight
Instrument Trays			
IT1015	9710015001	15" W x 10" D	3 lbs.
IT1621	9716002100	21" W x 16" D	4 lbs.
IT1525	3603915250	25" W x 16" D	5 lbs.
IT2025	3603920251	25" W x 20" D	6 lbs.
IT1219	9719126250	19" W x 12 1/2" D	4 lbs.
IT2520	9736392025	20" W x 25" D	6 lbs.
Catheter Tray			
CT-1	9785030101	8 1/2" W x 3" D x 1 1/2" H	1 lbs.
CT-IT	9782831000	8 1/2" W x 3" D x 1" H	1 lbs.
CT-IC	9782832000	8 1/2" W x 3" D x 1 1/2" H	1 lbs.



Iodine Cups
IC1-2



Medicine Cup
MC-1



Graduated
Measuring Cups
GM1-3

Model #	Item #	Capacity/Dimensions	Shipping Weight
Iodine Cups			
IC-1	9775500002	6 oz.	1 lbs.
IC-2	9710565001	14 oz.	1 lbs.
Medicine Cup			
MC-1	9784920001	2 oz.	1 lbs.
Graduated Measuring Cups			
GM-1	9795160001	16 oz.	1 lbs.
GM-2	9795320002	32 oz.	1 lbs.
GM-3	9795640003	64 oz.	1 lbs.



Surgical Utensils

Construction Specifications

Construction Specifications

Standard Dimensions	N/A
Material	304 stainless steel
Gauge	20GA: Graduated Measuring Cup, Instrument trays (exception 15"x10" Instrument Tray is 24GA) 22GA: Dressing Jars, Emesis Basins, Ointment Jar, Bucket Pail, Solution Basin, Catheter Tray and Cover, Iodine Cups, Graduated Measuring Cups
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	N/A

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications service for estimated shipping costs at info@blickman.com.

Toll Free: (800) 247-5070
Local: (973) 330-0557



www.blickman.com
info@blickman.com

1099 Wall Street West, Suite 250A, Lyndhurst, NJ 07071